



PURELINE OC8000 Oxygen Concentrator

SAVINGS FOR YOU, SAFETY FOR YOUR PATIENTS

Stop paying for tanks, refills, and deliveries.
Stop wasting time on tank billing, scheduling,
management, and maintenance.

INNOVATIVE CONCENTRATOR TECHNOLOGY

Pureline has perfected oxygen concentrator technology, added safety features like O2 flow and purity alarms.

The Pureline oxygen concentrator is UL tested and backed by the industry's best warranty.



CALCULATE SAVINGS
Visit www.superavet.com/ROI to
calculate YOUR actual savings!



Key Features:

- Extracts medical-grade oxygen (USP93%)
- Delivers on demand up to 8 lpm
- Connects with industry standard DISS fittings
- Includes 10' high-pressure green oxygen hose
- Quiet operation at 53 dBA
- Light weight 54 lbs.
- Energy efficient
- Built specifically for anesthesia machine use
- FDA approved, UL tested
- Meets 13485:2016 ISO & 60601-1 IE Standards
- Oxygen purity safety alarm
- Made in the USA

SUPERA

ANESTHESIA
INNOVATIONS

LEARN HOW PURELINE CAN BENEFIT ANY VETERINARY PRACTICE AT WWW.SUPERAVET.COM
(503) 723-5068 OR (877) 620-1500 INFO@SUPERAVET.COM

ANESTHESIA INNOVATIONS SINCE 1993

©2021 SUPERA ANESTHESIA INNOVATIONS. ALL RIGHTS RESERVED. PRODUCED IN USA.