

M5400 Mobile Anesthesia Cart

A convenient mobile platform for transporting or temporary anesthesia procedures. Large heavy duty 20" X 44" table top



Features:

- Heavy gauge medical grade stainless steel
- Raised drip edges to control fluids
- Large 5" locking casters
- Patented safety pop off valve - it saves patient lives!
- 0-4 Lpm oxygen flowmeter with mechanical stop control
- Oxygen flush with safety restrictor and button shield
- CO2 absorber with easy change "One Seal" canister
- Dual O2 supply ready with check valves
- Easy to use non-rebreathing circuit common outlet
- Syringe holder (for inflating E/T cuffs)

Dimensions: 20" X 44" X 36" high

Available in non-rebreathing P/N M5450

*Shown with optional vaporizer, e-tank manifold and waste gas interface.



SUPERA

ANESTHESIA
INNOVATIONS

Formerly LEI Medical
(503) 723-5068 or
(877) 620-1500 toll-free
info@superavet.com
12330 S.E. Highway 212
Clackamas, OR 97015