

# Model M1200 Rebreathing “A” Frame Stand Veterinary Anesthesia Machine

*Superior stability with our “A” Frame design stand*



*Note: Vaporizer & E-tank manifold are not included. Sold Separately*

## M1200 features:

- Patented Safety Pop-Off Valve
- Built in stainless steel I.V. pole
- 17 1/2” wide X 9” deep heavy duty monitor shelf
- 0-4 Lpm oxygen flowmeter with mechanical stop control
- Oxygen flush with safety restrictor and button shield
- Dual O2 supply input with check valves
- 15mm female common gas outlet for circuits
- Syringe holder (for inflating E/T cuffs)
- Mounting option for universal control arm for “Bain” circuits (see MA2013 below)



*I.V. Pole included*

*Optional Universal Control Arm*

*Model: MA2013*



## **SUPERA**

**ANESTHESIA  
INNOVATIONS**

**(503) 723-5068 or  
(877) 620-1500 toll-free**